

1

The Microfluidic probe: a non-contact scanning microfluidic technology

Emmanuel Delamarche, Robert Lovchik
Govind Kaigala, Julien Cors
<http://www.zurich.ibm.com/st/bioscience/>

© 2013 IBM Corporation

4

Principle of the Microfluidic Probe (MFP)

How does it work?

- The MFP is brought close to a surface in the presence of an immersion liquid
- a processing liquid is injected by one aperture and aspirated back together with some of the immersion liquid by a second aperture
- as a result, the processing liquid is focused giving the possibility of localizing (bio)chemical processes on a surface!

G. Kaigala, R. Lovchik, U. Drechsler, E. Delamarche, "A vertical microfluidic probe," *Langmuir*, 2, 2011, 5688-5693.
D. Juncker, H. Schmidt, E. Delamarche, "Multifunctional microfluidic probe," *Nature Materials*, 4, 2005, 622-628.

© 2013 IBM Corporation

7

Our specific MFP implementation

Key points

- compatible to standard silicon manufacturing processes
- 100s of MFP heads / wafer
- many designs co-fabricated
- reasonably easy to set-up and use system
- requires standard inverted microscope

G. Kaigala, R. Lovchik, U. Drechsler, E. Delamarche, "A vertical microfluidic probe," *Langmuir*, 2, 2011, 5688-5693.

© 2013 IBM Corporation

2

Experimental Biosciences @ IBM Research – Zurich

Chips for Life Brain Chip Microfluidic probe

IBM Research - Zurich Nanotech Center

© 2013 IBM Corporation

5

Video showing the confinement as a function of Qa/Qi

Hydrodynamic Flow Confinement in a 10-µm-Deep Gap

Real time video taken through a glass substrate using an inverted microscope

video available at:
<http://www.zurich.ibm.com/st/video/mmat1435-s1.mov>

© 2013 IBM Corporation

8

Microfabrication of vertical MFP heads

- spin coat resist
- expose and develop
- CMP
- remove resist
- spin glass wafer
- bond wafer to Si
- Etching of channels with mask at water level
- Etch and separate MFP
- spin and polish Si wafer
- remove any wax with solvent

G. Kaigala, R. Lovchik, U. Drechsler, E. Delamarche, "The vertical microfluidic probe," *Langmuir*, 2011, 2, 5686-5693.

© 2013 IBM Corporation

3

What is special about the microfluidic probe?

"classical" closed microfluidic system The Microfluidic Probe (MFP): a non-contact, scanning microfluidic system

Web of Science® Results: 4 (including this article) Multiple microfluidic probes

Stephen R Quake et al. Nature Biotechnology 22, 435-439 (2004) D. Juncker, H. Schmidt, E. Delamarche, Nature Materials, 2005, 4, 622-628.

© 2013 IBM Corporation

6

MFP – Research Landscape

Single devices multiple apertures devices with windows vertical MFP

IBM Research

Olav Öster and Gösta Johansson, Chalmers, Sweden (2) microfluidic pipette Tomokazu Matsuda, Tohoku University, Japan (1)

1. H. Shiku et al., "Microfluidic dial-facility probe to collect messenger RNA from adherent cells and spheroids," *Analytical Biochemistry*, 2009, 385, 135-142.
2. A. Kishi et al., "A Microfluidic Pipette for Single-Cell Pharmacology," *Analytical Chemistry*, 2010, 82, 4523-4530.
3. M. A. Gesswein, T. Gervais, and D. Juncker, "Microfluidic quadrupole and floating concentration gradient," *Nature Commun.*, 2, 2011, 454.

© 2013 IBM Corporation

9

Protecting channels during fabrication – key challenge!

(a) channel protection with wax

(b) Wax removal

- For improved yield in fabrication wax is filled in the channels to minimize clogging of with debris generated during fabrication and packaging (dicing, polishing)
- Wax could be left in the channel post-processing for packaging and shipping of the MFP heads.

© 2013 IBM Corporation

10 MFP head - Advanced designs

G. Katgika, R. Lovchik, U. Drelich, E. Ddamache, "The vertical microfluidic probe," *Largomir*, 2011, 2, 5686-5693.

12 Microfluidic quadrupole and floating concentration gradient

- Possibility to work with various environments
- Spatial multiplexing of liquids
- Design for a 'local minima' within the flows for trapping particles

M. A. Qasimab and D. Jenckel, 1844-1846, *Proceedings microTAS*, 2010
M. A. Qasimab, T. Garvan, D. Jenckel, *Nat. Commun* 2011, 2, 464.

16 Patterning and removing Proteins onto from a Surface

- surface-density gradients of proteins on a glass surface
- gradient was formed by varying the writing speed of the MFP
- fluorescently labeled proteins removed from a surface
- the processing liquid contained a surfactant, a high pH and high ionic strength
- proteins adsorbed on a glass slide are removed by the processing liquid

11 Parameters determining the hydrodynamic flow confinement

variable parameters:

- Q_1 and Q_2 (nL/min)
- gap between the MFP and substrate (μm)
- writing speed (mm/s)

apertures: 20 μm
distance between the apertures: 120 μm
gap: 2 μm

14 Writing protein patterns using a MFP

- Proteins in the processing liquid deposit on the scanned surface
- no drying artifact due to the presence of the immersion liquid (biological buffer)
- moving quickly the MFP makes the immersion liquid inserting below the processing liquid → **non-writing mode!**

17 Staining cells & picking cells

12 Parameters determining the hydrodynamic flow confinement

- fixed parameters
 - size and shape of the apertures
 - number of apertures
 - relative position of the apertures

MFP with 2 concentric apertures

15 Patterning proteins on a surface using a MFP

Fluorescence microscope image

- 2 types of antibodies were subsequently patterned on an activated glass slide
- array has 1384 spots spaced 80 μm apart
- ~130 μL of antibody solution and 0.3 s writing time per spot
- array needed 300 nL of antibody solution and 15 min writing time
- deposited, removing, delivering and creating protein gradients on surfaces in a way that preserves the protein function
- functionalize surfaces to study advanced cell-surface interactions

18 Microscale tissue staining

19

Tissue analysis: Needs and challenges

Tissue analysis is critical for proper tumor identification and therapies!

Strong desire for more and better "quality" information

- multiplexing
- adaptive staining conditions
- use/reed for lesser amount of tissue
- supply data for complex decision algorithms
- use/reed for lesser amount of reagents (antibodies)

Need for more sample

H. Ahmadzadeh, L. V. Thompson and E. A. Arlt, *Anal Bioanal Chem*, 2006, 389.

© 2013 IBM Corporation

22

Flow stability and staining contrast

- flow confinement is stable for tens of minutes with area variations < 5% of its mean
- tissue section is not delaminated from glass slide by shear stress from the flow

Human thyroid (normal) tissue section, 4 μm thick, processed with anti-thyroglobulin (anti-TGB)

© 2013 IBM Corporation

25

Tissue Microarrays

- breast cancer tissue cores (30 individuals)

Core	Stage	Grade	ER	PR	p53
A1	IA	0	+	+	+
B1	IA	0	+	+	+
B2	IB	1	+	+	+
C1	IIA	1	+	+	+

R. D. Lovchik, G. V. Kagalja, M. Georgiadis, E. Delamaro, "Micro-immunohistochemistry with the microfluidic probe" Lab on a Chip, 12, 2012, 1040-1043.

© 2013 IBM Corporation

20

Miniaturized Immunohistochemistry: State-of-the-art

Microfluidic immunohistochemistry platform

PARTECELL: partial treatment of cells using laminar flows

Layered immunohistochemistry

A. T. Chin, M. A. M. Cui et al. Proceedings of AP International Conference 2012, 1-11, 1988-2012.

G. M. Whiteside, J. Clin. Invest. 118, 2008.

Bo Song, Manli A. M. Gloger and Jean-Charles G. Stroh et al., Analyst, 134(10) 2005, 2006.

© 2013 IBM Corporation

23

Flow stability on a human cancerous thyroid tissue section

© 2013 IBM Corporation

26

Tumor tissue processing using the vMFP

- breast cancer tissue sections – anti-p53 antibodies (spots), anti-progesterone receptor (PR) antibodies and hematoxylin line

R. D. Lovchik, G. V. Kagalja, M. Georgiadis, E. Delamaro, "Micro-immunohistochemistry with the microfluidic probe" Lab on a Chip, 12, 2012, 1040-1043.

© 2013 IBM Corporation

21

First steps with local staining: Counter staining using a MFP

stain	target
DAPI	nucleus
phalloidin	actin
hematoxylin	nucleus
eosin	nucleus
	cytoplasm

breast tissue

100 μm

© 2013 IBM Corporation

24

Micro-Immunohistochemistry (staining)

1. TISSUE SELECTION
2. DEWAXING & REHYDRATION
3. LOCAL PROCESSING
4. POST-PROCESSING

reduced processing/incubation time probably due to increased connection

R. D. Lovchik, G. V. Kagalja, M. Georgiadis, E. Delamaro, "Micro-immunohistochemistry with the microfluidic probe" Lab on a Chip, 12, 2012, 1040-1043.

© 2013 IBM Corporation

27



Local dewaxing

- key goal: local staining AND local dewaxing
- a sample is not "spent" after a single analysis
- enabling retrospective analysis taking advantage of the >1.5 billion tissue samples in bio-banks?

MFP head designs and local dewaxing

G. Kagalja et al., pTAS 2012 | October 29, 2012

© 2013 IBM Corporation



30

Summary

- The MFP is a non-contact microfluidic technology; it works in the presence of a surrounding liquid (no drying artifacts)
- it can write patterns with $\sim 5 \mu\text{m}$ accuracy with additive and subtractive surface patterning processes
- it is very versatile
- it represents an exciting opportunity for microscale tissue section processing

Acknowledgements

Marios Georgiadis (now at ETH-Zurich), Martina Hitzbeck, Ute Drechsler, Michel Despont and Walter Riess

IBM (\$\$) and "BioProbe" ERC grant (EU) (\$)

© 2013 IBM Corporation